



Testosterone Metabolites Profile

COLLECTION & SHIPPING INSTRUCTIONS

Androstenedione
Testosterone
5 α -dihydrotestosterone (DHT)
5 α -androstan-3 β ,17 β -diol (3 β -adiol)
5 α -androstan-3 α ,17 β -diol (3 α -adiol)

KIT CONTENTS

- | | |
|---------------------------|--|
| 1 - Requisition Form | 1 - Ice Pack |
| 1 - Instruction Sheet | 1 - Biohazard Bag |
| 1 - Questionnaire | 1 - Styrofoam Container |
| 2 - Gel-Free Red Top Tube | 1 - Cardboard Shipping Box |
| 1 - Transfer Tube | 1 - UPS Next Day Air Pack (US & Canada only) |
| 1 - Disposable Pipette | 1 - UPS Next Day Shipping Label (US & Canada only) |

CAREFULLY READ COLLECTION AND SHIPPING INSTRUCTIONS THOROUGHLY BEFORE STARTING.
UPON RECEIPT – REMOVE LID FROM STYROFOAM CONTAINER AND PLACE INTO THE FREEZER – **DO NOT REMOVE THE ICE PACK.**

PLEASE NOTE: When shipping the specimens back to the lab, send Monday through Thursday Only. **DO NOT SHIP ON FRIDAYS AS NEXT DAY DELIVERIES ARE NOT ACCEPTED ON SATURDAY OR SUNDAY.**

OBSERVED HOLIDAYS ARE AS FOLLOWS: NEW YEARS DAY, MEMORIAL DAY, FOURTH OF JULY, LABOR DAY, THANKSGIVING DAY AND THE DAY AFTER, AND CHRISTMAS DAY.
DO NOT SEND SAMPLE(S) THE DAY PRIOR TO OR THE DAY OF A HOLIDAY.

COLLECTION INSTRUCTIONS

- *PATIENT MUST BE FASTING FOR 8-12 HOURS PRIOR COLLECTION**
- *BLOOD DRAW MUST OCCUR BEFORE 10AM**
- **HEMOLYSIS OR RED BLOOD CELLS IN SERUM WILL BE REJECTED****

1. Draw a **gel-free** red top tube completely; allow at least 15- 30 minutes for the blood to clot.
2. Centrifuge for at least 10 minutes.
3. Pipette the separated serum into the transfer tube. We need a minimum of 4mL of serum. ***Do not pick up any Red blood cells.***

Please call 206.209.4200 with questions or Toll Free 855.405.8378
Monday through Friday 6:00 am – 6:00 pm PST

6839 Fort Dent Way, Ste 206
Tukwila, WA 98188
Fax: 206.209.4211
www.meridianvalleylab.com info@meridianvalleylab.com

4. Label the transfer tube as **SERUM** with patient's name and date collected.
5. Freeze the serum tube immediately. Keep frozen until ready to ship.
6. Complete the enclosed requisition form; include payment (if required)

SHIPPING INSTRUCTIONS

1. Put the requisition form into the pocket on the outside of the biohazard bag.
2. Place the frozen transfer tube in the biohazard bag and seal.
3. Place the biohazard bag next to the frozen ice pack in the Styrofoam container.
4. Place the Styrofoam container in the cardboard shipping box.
5. Place the shipping box in the UPS next day air pack and seal securely.
6. Place return shipping label with barcode on the outside UPS next day air bag.
7. Save return receipt with the tracking number for your records.
8. Send specimens Monday through Thursday only.
9. Call UPS toll free at **1.800.742.5877** for a pickup or for the nearest drop location. UPS will only pick up on weekdays.
10. When calling for UPS pickup, request "ON CALL AIR" only. **DO NOT request a one-time pickup.** Meridian Valley Lab pays for "ON CALL AIR" only. UPS expects the patient to pay for a one-time pickup.

Please call 206.209.4200 with questions or Toll Free 855.405.8378
Monday through Friday 6:00 am – 6:00 pm PST

6839 Fort Dent Way, Ste 206
Tukwila, WA 98188
Fax: 206.209.4211
www.meridianvalleylab.com info@meridianvalleylab.com

Male Serum Testosterone Metabolites Profile (4416, 4417)

Hormone Symptoms Questionnaires

Last Name: _____ First Name: _____

Age: _____ Height: _____ Weight: _____ Blood Pressure: _____

Previous / Current use of Hormones / Medications: Please indicate any hormone(s) you have used in the past **2 months** as shown in the example below. Also list other medications or herbal supplements you are taking. (Use back if necessary)

FOR LAB USE ONLY. THIS WILL ONLY BE USED INTERNALLY AND WILL BE DISCUSSED ONLY WITH YOU OR YOUR DOCTOR.

Hormone Therapies	Testosterone	DHEA	5a-reductase inhibitor	hCG	
Hormone Type					
Brand Used					
Delivery					
Date- last used					
Time -last used					
# of times/day					
How long used					
Other Medications					

Symptom: Please indicate the symptoms you are experiencing as **0 (none)**, **1 (mild)**, **2 (moderate)**, **3 (severe)** by darkening the number in the box.

<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Burned Out Feeling	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Muscle Size	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Cold Body Temperature
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Mental Sharpness	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Neck Or Back Pain	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Libido
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Increased Forgetfulness	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Sore Muscles	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Erections
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Apathy	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Stamina	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Prostate Problems
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Depressed	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Flexibility	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Increased Urinary Urge
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Nervousness	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Increased Joint Pain	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Decreased Urine Flow
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Anxious	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Bone Loss	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Male Pattern Baldness
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Irritable	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Heart Diseases	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	ringing In Ears
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Chronic Stress	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Heart Palpitations	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Dizzy Spells
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Acute Stress	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Elevated Cholesterol	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Allergies
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Difficulty Sleeping	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Elevated Triglycerides	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Sensitivity To Chemicals
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Morning Fatigue	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Diabetes	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Headaches
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Evening Fatigue	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Weight Gain - Waist	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Other _____
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Night Sweats	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Weight Gain - Breast	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	
<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Hot Flashes	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	Sugar Craving	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3	