

6839 Fort Dent Way, Ste. 206 Tukwila, WA 98188 tel 206.209.4200 • 855.405.TEST (8378) Fax 206.209.4211

Instructions For	Dried Urine Hormone Test(s)
Notes for this test:	Do not consume creatine or supplements containing creatine during the collection period. If you are using estrogen (bi-est, estradiol, estriol, etc), testosterone or cortisol cream vaginally or in the labia region, please indicate that on the requisition form and the hormone questionnaire.

Kit Contents					
1 Requisition Form	1 Questionnaire				
1 Instruction Sheet	1 Collection Cup				
4 Dry Urine Cards	1 Biohazard Bag				
4 Green Stickers	1 Return Envelope				

Collection Notes

- Fluid intake volume and pattern can affect hormone concentration in urine. In order to maintain a relatively steady urine concentration, avoid consuming large volumes of fluid all at one time. Instead, sip throughout the day. Totally fluid intake on during the collection period should not exceed 3 liters.
- Avoid supplements containing creatine the day before and during the collection period.
- If on hormone replacement therapy, (estrogen, testosterone, thyroid, etc.) *consult with your physician* as to whether these hormones should be taken during collection. It is usually recommended to continue hormone replacement during collection.
- Women with regular menstrual cycles should collect urine on day 19, 20, or 21 of a 28-day menstrual cycle.
 (Day 1 of the menstrual cycle is the first day of bleeding.) The monthly cycle is longer or shorter than 28 days, add or subtract a corresponding number of days and adjust the collection date.
- Postmenopausal women using estrogen and progesterone should collect when both estrogen and progesterone have been used for 5 continuous days. Postmenopausal women not using hormones may collect urine on any day.

Before Collection -

Before starting collection, review all of the collection instructions and timing instructions. Try to plan your day, as much as possible, based on the collection timing.

We recommend using an alarm, cell phone, watch, or other timing device to help alert you of collection times to ensure you do not miss one.

Hormone Replacement Users Please Note

If you are using estrogen (bi-est, estradiol, estriol, etc), testosterone or cortisol cream vaginally or in the labia region, please indicate that on the requisition form and the hormone questionnaire.



Fill out all information on the cards. Collect urine into cup. Make sure information is filled out prior to starting.



Dip the card, submerge once or twice, to just above the top marker line. Discard urine after each collection



After removing the card from the urine cup, tape it so that it hangs freely to dry for at least 24 hours.



Repeat for each urine collection according to the collection timing guidelines.

Collection Timing

Card #1 First morning urine (5:00-8:00 am): Collect immediately upon waking. After collection, you may take supplements, medication, and have breakfast as usual. After first morning collection, limit total fluid intake to 8-10 oz until second collection.

Card #2 Second morning urine collection (7:00-11:00 am): Collect at least 2 hours after previous collection.

Card #3 Afternoon collection, before dinner (3:00 -6:00 pm): Collect at least 4 hours after previous collection. Restrict fluid intake 2 hours before collection. If you cannot go without fluid for 2 hours, take small sips of water.

Card #4 Bedtime collection (9:00 pm— Midnight): Collect at least 4 hours after previous collection. Restrict fluid intake 2 hours before collection.

Allow cards to dry in a well-ventilated area. Use the provided green sticker to hang them from a counter edge or towel bar. Do not allow the filter paper portion of the card to come into contact with anything while drying. Dry the cards for at least 24 hours.

Once samples are dry, fold outer cover back over urine strip and place into plastic bag, taking care not to touch the collection paper portion.

Complete the requisition form along with payment (if necessary).

Remember to fill out the enclosed questionnaire.

Shipping Instructions

Place the cards into the biohazard bag.

Fold the requisition form and place it in outside pocket of the biohazard bag.

Place the biohazard bag into the provided return envelope.

Contact the provided courier service, or drop-off at an approved drop zone/facility.



0 1 2 3 Aches And Pains

0 1 2 3 Loss Scalp Hair

0 1 2 3 Tender Breasts

0 1 2 3 Sensitivity To Chemicals

0 1 2 3 Night Sweats

0 1 2 3 Memory Lapse

0 1 2 3 Headaches

0 1 2 3 Fibromyalgia

0 1 2 3 Sugar Craving

0 1 2 3 Stress

0 1 2 3 Loss of Eyebrows

0 1 2 3 Bleeding Changes

0 1 2 3 Increased Facial And Body Hair

0 1 2 3 Elevated Triglycerides

6839 Fo	rt Dent Way, Suite 206
	Tukwila, WA 98188
206.209.4200 •	855.405.TEST (8378)
	fay 206 209 4211

	• • • • • • • • • • • • • • • • • • • •	ey LAB	200	stionna			
Last Nar	ne:			_ First	Name:		
Age:	Height:	W	eight:	Blood Pressure:			
the examp	ole below. Also list o	f Hormones / Medic ther medications or he WILL ONLY BE USED	rbal supplements yo	u are taking. (U	lse back if neces	sary)	t 2 months as shown in
Prescriptive Meds			1	2		3	4
Medication	n	Progesterone	Estrogen	Testos	sterone	Cortisol	
Brand Use	ed						
Delivery (oral, cream, subling	ual)					
Amount (i	n mg's)						
Date and	time last used						
# of times	/days, days/month						
How long	used						
Other Sup	pplements	(Example)					
	•	DHEA/Pregnenolone					
		Soy / Flaxseed					
		DIM / I3C					
Please	Report the sympt	oms you are current	Male Que	estionna	ire		
	Burned Out Feeling	0 1 2 3 Apathy		0 1 2 3 Difficu	ulty Sleeping	0 1 2 3	Increased Forgetfulness
	•	arpness 0 1 2 3 Depres			al Fatigue		Irritable
	Nervousness	0 1 2 3 Anxiou			ng Fatigue		Evening Fatigue
	Decreased Stamina		sed Muscle Size		Muscles	0 1 2 3	Increased Joint Pain
0 1 2 3	Decreased Flexibility		r Back Pain	0 1 2 3 Weigl	ht Gain - Breast or I	Hips 0 1 2 3	Weight Gain - Waist
0 1 2 3	Elevated Triglycerides	0 1 2 3 Sugar	Craving	0 1 2 3 Heart	Palpitations	0 1 2 3	Dizzy Spells
0 1 2 3	Headaches	0 1 2 3 Ringing	ı In Ears	0 1 2 3 Cold I	Body Tempurature	0 1 2 3	Allergies
0 1 2 3	Sensivitivity To Chemi	cals 0 1 2 3 Decrea	sed Erections	0 1 2 3 Decre	eased Libido	0 1 2 3	Prostate Problems
0 1 2 3	Decreased Urine Flow	0 1 2 3 Increas	ed Urinary Urge	0 1 2 3 Stres	s	0 1 2 3	Other
			Female Q	uestionn	aire		
	Regular Cycles	Hysterecto	omy:	☐ No ☐ Ye	es Ova	aries Removed:	☐ No ☐ On ☐ Both
j	Irregular Cycles	Pregnant:		□ No □ Ye	es Month of Preg	nancy	
ļ	☐ No menstrual Cycle	es Polycystic	Ovarian Syndrome	□ No □ Ye	es		
Date of Las	t Menses	For specimen of	ollection: # of days in c	ycle from day 1 of	last menses		
When was	the last time you used	hormone based birth contro	ol (pills, IUD, etc.)	H	low long were you o	on it?	
0 1 2 3	Hot Flashes	0 1 2 3 Bone L	oss	0 1 2 3 Wate	r Retention	0 1 2 3	Acne
0 1 2 3	Vaginal Dryness	0 1 2 3 Depres	sed	0 1 2 3 Tearf	ul	0 1 2 3	Irritable
0 1 2 3	Foggy Thinking	0 1 2 3 Anxiou	s	0 1 2 3 Morni	ng Fatigue	0 1 2 3	Evening Fatigue

0 1 2 3 Sleep Disturbed

0 1 2 3 Weight Gain - Hips

0 1 2 3 Heart Palpitations

0 1 2 3 Decreased Libido

0 1 2 3 Incontinence

0 1 2 3 Cold Body Temperature

0 1 2 3 Nervousness

0 1 2 3 Mood Swings

0 1 2 3 Uterine Fibroids

0 1 2 3 Fibrocystic Breasts

0 1 2 3 Allergies

0 1 2 3 Other _

0 1 2 3 Weight Gain - Waist